SprintRay EU Denture Base Instruction for use

1. Indication

SprintRay EU Denture Base is a light-curable, polymerizable resin to be used in conjunction with extraoral light-curing equipment.

SprintRay EU Denture Base is indicated for the fabrication and repair of removable full and partial dentures and base plates by additive manufacturing.

2. Contraindication

SprintRay EU Denture Base

restorations are contraindicated ...

1. ... if a patient is known to be allergic to any of the ingredients.
2. ... direct intraoral contact of uncured or partially uncured material.
3. ... for every application that is not part of the indication (see above).

3. Requirements

Software:

- RayWare (CAD)

Hardware:

- Pro Desktop 3D (3D Printing)

- Pro Wash / Dry (Cleaning)

- ProCure (Post Curing)

SprintRay, Inc., 3577 North Figueroa St., Los Angeles, CA 90065

4. Material

SprintRay EU Denture Base consists of functional (meth) acrylic resins, initiators, dyes and stabilisers.

5. Geometric presetting

Lingual Ridge ≥ 2.5 mm (Maxilla / Mandibula)

Palatal/Lingual ≥ 3.0 mm (Maxilla) / ≥ 2.5 (Mandibula)

Facial/Buccal ≥ 2.5 mm (Maxilla / Mandibula)

6. Material parameters

Radiation penetration depth controllable by exposure time

≥ 50 µm (4 mils)

≥ 100 µm (5 mils)

7. Manufacturing process (fig. 1-10)

1. Prepare data (RayWare).
2. Choose process parameters (build style, etc.).
3. Transfer prepared data to 3D printing device.
4. Prepare 3D printing – shake bottle.
5. Fill resin tank of 3D printing device.
6. Build the parts.
7. Clean parts (ProWash / Dry): 5 min. / 5 min.
8. Dry parts (ProWash / Dry): 3 min.
9. Post curing (ProCure): 5 min. / 60 °C.
10. Finish parts.

8. Bonding of denture teeth

If using printed denture teeth – see additional IFU **SprintRay EU Temporary Crown & Tooth** for details. In order to realize the optimal bonding of artificial teeth to the dental object, it can be necessary to roughen the surface of the denture and to use an appropriate adhesive. The suitability of a bonding agent or adhesive must be approved prior use. Recommended bonding agent is VITA VIONIC® BOND. The commonly used techniques and dental work steps must be carefully considered when manufacturing the dentures. Please follow the instructions and indications of the corresponding suppliers.

9. Notice

Please follow the instructions provided by the software manufacturer in terms of parameter settings and design recommendations. Please follow the instructions provided by the hardware manufacturer in terms of parameter settings/printing & post curing recommendations. To avoid detrimental effects on material quality do not expose the liquid material to irradiation under any circumstances. Deviations from the described manufacturing processes or storing conditions may lead to different mechanical and optical properties of the material. Ensure personal protective gear during processing.


Caution: Polymerized resins are chemically resistant - avoid stains on clothing! Avoid any contact with skin and eyes! In case of accidental contact, rinse with adequate running water, consult a doctor if necessary. The lot number and the best before date are indicated on each packaging. In case of claims please always indicate the lot number of the product. Do not use the product after expiry of the best before date. Dispose of contents/container according to official regulations.

10. Hazard statements

Causes skin irritation. May cause an allergic skin reaction.

Causes serious eye irritation. May cause respiratory irritation.

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light-pink/
hell-pink/
rose-clair/
rosa-claro/
rosa-chiaro
REF: SRE-1001303

deep-pink/
dunkel-pink/
rose-profond/
rosa-profundo/
rosa-scuro
REF: SRE-1001304

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