

0 mm Light tested directly on sensor

Light Mode: Standard

Irradiance Measurement at 0 mm

1321* mW/cm²

6 mm Light tested 6 mm from sensor

Light Mode: Standard

Irradiance Measurement at 6 mm

1166* mW/cm²

6 mm distance using tooth



Elipar DeepCure-L

3M

Light Type	LED-A
Stated Irradiance	1470 mW/cm ²
Variance	-10%
Light Tip	Light Guide 10 mm
Modes Available	1

Notes

Damage	<input type="checkbox"/>
Debris	<input type="checkbox"/>
Barrier	<input type="checkbox"/>
Demo	<input type="checkbox"/>

Manufacturer/Brand

Category

Filtek Bulk Fill Flowable / 3M

U

Filtek Bulk Fill Flowable / 3M

A1, A2, A3

Filtek Supreme Ultra Flowable / 3M

Opaque shades

Filtek Supreme Ultra Flowable / 3M

All other shades

Min Curing Times

Max Increment

0 mm

6 mm

4.0

10

10

4.0

20

20

1.5

20

20

2.0

10

10

Notes

Potential Heat Concern

14 s *

Distance and Area can have a significant effect on irradiance.