

3Shape TRIOS

Better for them.
Better for you.



3shape 

Welcome to the world of TRIOS digital dentistry.

Dental treatments that reflect your standards.

You can rely on every dental impression taken with a TRIOS® being at least as accurate as a conventional impression*.

Learn more at 3Shape.com/treatments

Patient engagement you can sense.

Boost case acceptance by engaging with your patients. Together with the 3D color scans and the TRIOS engagement apps, they'll notice the difference.

Learn more at 3Shape.com/engagement

An experience that just feels right.

Your digital experience should make you smile. So we've partnered with 25,000+ labs and providers to make it easy for you to find the one you want.

Learn more at 3Shape.com/experience

A digital experience you can be proud of.

Your patients will be happier and your practice more successful with a TRIOS, and we'll make sure you have everything you need to get started.

Learn more at 3Shape.com/practice





3Shape TRIOS

The world's most recognized intraoral scanner brand* – perfected over 5 generations.

That's five generations of pushing the limits of innovation in digital dentistry. Of surpassing expectations to deliver ease-of-use, award-winning design, and peace of mind.

And of breaking barriers by seamlessly connecting dentists with leading solutions and labs through the dental world's most collaborative platform: 3Shape Unite.

Learn more at 3Shape.com/trios

Unparalleled dental ecosystem.

The 3Shape Unite workflow engine unlocks a world of possibilities with your scan. Utilize thousands of dental apps, simplified treatment workflows, and access to your TRIOS scans in the cloud from any device.



Accuracy with realistic color.

With any TRIOS scanner, it's a breeze to create accurate scans in realistic colors and shade measurements. Impressing not only your patients but making it easier than ever for your lab to match colors.



Dynamic occlusion.

We've found a way for you to provide even more information in your restorative cases. How? With TRIOS Patient Specific Motion. Capture the bite with more occlusion information, which may improve the fit of your restorations and make your patient even happier with their outcomes.



Which TRIOS scanner is right for you?

There's a 3Shape TRIOS for everyone. Whichever model you choose, you benefit from industry-leading scanning experiences, and unparalleled clinic-to-lab workflows.



TRIOS 5 Wireless Intelligent. The ultimate TRIOS.

Next-level ergonomics.
Effortless scan technology.
Hygienic by design.

TRIOS 4 Wireless Advanced. Extra potential.

2-in-1 wireless technology.
Detect wear and tear
Scan-ready in seconds.

TRIOS 3 Wired Proven. An industry winner

Trusted technology.
Invisalign connected.
Improved experience

TRIOS Core Wired Core. Essential functionality.

Core scan technology
Built for send-to-lab

Proven clinical accuracy	✓	✓	✓	✓
Realistic color scans	✓	✓	✓	✓
Automatic scan clean up	✓	✓	✓	✓
Access to Unite ecosystem	✓	✓	✓	✓
Print-ready model design	✓	✓	✓	✓
Patient engagement apps	✓	✓	✓	
MOVE+ compatibility	✓	✓	✓	
Full-mobility wireless scanning	✓	2-in-1 SmartCord		
TRIOS Share capability	✓	✓		
Preventative care with surface caries detection	✓	✓		
Enhanced tip features	Calibration-free scanning Option of TRIOS Ready Tip	Smart tip calibration counter InstantHeat technology		
Smart battery charging	✓			
Energy-saving sleep mode	✓			
ScanAssist intelligent-alignment technology	✓			
Sensory feedback guides scanning	✓			
Hygienic design	Sleeve and TRIOS Ready Tip			

TRIOS 5. The ultimate TRIOS.

Intraoral scanning designed for an effortless experience, next-level ergonomics, and improved hygiene.

Learn more at [3Shape.com/trios5](https://3shape.com/trios5)



Ergonomic.

30% smaller and 20% lighter than TRIOS 4. Our smallest and lightest scanner to date. Made to fit perfectly in your hand.

Effortless.


ScanAssist makes TRIOS 5 our easiest-to-use. LED and haptic feedback guide you as you scan. One battery charge will last you all day.

Hygienic.

Re-designed autoclavable scan tips and body sleeves significantly reduce risk of patient cross-contamination. Hygiene by design.

Convenient.

Discover convenience redefined with a single-use TRIOS Ready Tip that's always ready to use, straight from the box. With the Ready Tips, your team can spend less time on cleaning, and more time on patient care.



TRIOS 4. Advanced potential.

Advanced wireless scan technology expanded for even better treatment options.

Learn more at 3Shape.com/trios4



Wired and wireless in one.

Two-in-one scanning combines the freedom of wireless dentistry with the reassurance of a cord. All with the convenience of a single scanner.

Detect wear and tear.

TRIOS Patient Monitoring helps you compare and track patient cases from previous visits, so you and your patient can see what's happening on screen.

Scan ready in seconds.

With the use of InstantHeat technology, our 4th generation scanner tips ensure the mirror inside the tip heats in a brief moment and no fog will compromise your scan accuracy (in known conditions).

TRIOS 3. The trusted choice.


Proven scanning technology for an even greater patient experience.

Learn more at 3Shape.com/trios3



Improved experience.

30% lighter. 10% smaller. The improved TRIOS 3 pod takes up less of your precious workspace. It has fewer cables, for a lot less fuss. And the LED indicator light lets you know when it's connected and ready to go!




New ways to engage patients.

TRIOS 3 also comes with free engagement apps. There are so many new ways to engage with your patients, add weight to your treatment proposals and gain – even better – case acceptance.

Trusted scanning.

Tens of thousands of dentists use – and love – award-winning TRIOS 3. It's a reliable, AI-driven tool that makes it a breeze to create accurate scans, in realistic colors, and send them off to the lab of your choice. Just like that!



Welcome TRIOS Core! Essential scan technology.

Essential scan technology built for send-to-lab,
with TRIOS accuracy as standard.

Learn more at 3Shape.com/trios-core



Never compromise on the essentials.

With TRIOS Core, you get TRIOS scan technology in its most affordable form. AI-powered, TRIOS Core includes the essential scan functionality for accurate scanning, ease-of-use, and speed.




Built for send-to-lab.

The simple send-to-lab workflow is built to be easy and intuitive. Simply connect with your favorite lab, scan, and send. No file exports. No email attachments. Just a seamless time-saving workflow.

3Shape Unite. The TRIOS ecosystem.


A great scan is just the beginning,
unlock a world of possibilities.

Learn more at 3Shape.com/unite



Patient cases made available - from any device!

With the Unite Cloud app and web, store and manage patient scans, track case progress, communicate with your lab and view all your 3Shape app licenses directly in the cloud, no matter where you are.




Treatments made simple.

Experience guided and easy workflows to manage cases and communicate with your lab or preferred partner. Connectivity to more than 25,000 labs in just one click.



Dental apps made accessible.

Boost your work and collaborations with 3Shape, labs, and more than 100 world-class partners—from clear aligners, 3D printing and milling, patient management, and more. Use dental apps for routine, in-demand, and specialist treatments, all in one place.



3Shape Unite. Apps to treat and engage.



TRIOS Engagement Apps.

Our engagement apps make it easier than ever to communicate with your patients and help them understand their treatment needs.

Learn more at 3Shape.com/engagement-apps

TRIOS Patient Monitoring.

Accurately measure changes in dentition over time by scanning every patient at every appointment.

TRIOS Smile Design.

TRIOS Smile Design lets you digitally design a patient's potential new smile before their eyes.

TRIOS Treatment Simulator.

TRIOS Treatment Simulator simulates realistic corrections on-screen based on a 3Shape TRIOS intraoral scan.

TRIOS Studio Apps.

Together with our selection of Studio Apps, 3Shape TRIOS gives you the option of offering your patients fast, hassle-free, same-day treatment.

Learn more at 3Shape.com/studio-apps

TRIOS Implant Studio.

Plan, design and print custom surgical guides and temporary screw-retained crowns in one workflow.

TRIOS Design Studio.

Design and produce a broad range of restorations, including standard and screw-retained crowns.

TRIOS Splint Studio.

Design and produce splints, night guards and mouth guards with a guided, fast and intuitive workflow.

TRIOS Share.

Pass the scanner, share the power.

TRIOS Share technology enables you to digitize your entire clinic with just one TRIOS wireless intraoral scanner. Walk around with only one wireless TRIOS and use it on every PC in your practice via your Wi-Fi network.

Learn more at 3Shape.com/share



TRIOS MOVE+

It just got easier to engage.

With the Red Dot celebrated ergonomic design and compatibility with TRIOS 3, TRIOS 4 and TRIOS 5, MOVE+ is the perfect addition to your clinic. It's also the easiest way to move 3D digital scans and TRIOS engagement apps right where your patients can see them best - in the chair! TRIOS MOVE+ is easy to adjust, easy to use, and easy to move.

Learn more at 3Shape.com/triosmove




TRIOS Services. We've got you covered!

Start your TRIOS journey with ease! With personalized learning paths, you can progress at a pace that is right for you. And if you ever have questions, our support team is just a call or click away.

Plus, your TRIOS investment is protected by our top-notch service plan. 3Shape TRIOS Care means all-around coverage, less stress, and more focus on providing a great patient experience from day one.

Dive into tailored training, find answers whenever you need them, and enjoy the peace of mind that comes with our TRIOS Care service plan - first year is on us!

Learn more at 3Shape.com/trios-services



	With TRIOS you have access to	Add on top with TRIOS Care
TRAINING		
Get started with a tailored onboarding program and expert trainings		✓
Grow your skills with e-learning courses, webinars, or in-person training	✓	✓
Boost your business with insights into efficient practice management and growth	✓	✓
SUPPORT		
Unlimited support by 3Shape experts via phone, chat, or email		✓
Connect with 20,000+ dentists worldwide in our community forum	✓	✓
Troubleshoot from our library of self-guide help articles	✓	✓
INSURANCE		
5-year full warranty		✓
Accidental drop coverage with express replacement		✓
10TB case data cloud storage		✓

About 3Shape

3Shape was founded in 2000 by two people determined to unlock opportunities in 3D digital technology. From two, we've grown to over 2,000 – and from our beginning in Denmark, we're now present in over 100 countries around the world.

Our singular purpose remains the same – to advance and connect medical professionals and patients.

We drive progress in digital dentistry with groundbreaking, industry-recognized scanners for clinics and labs. And connect dentists and technicians with software workflows that just flow. We advance and coach professionals on every step of their digital journey – helping to truly elevate the customer experience.

Our impact is reflected in the results we achieve. Every second of every day, our award-winning products restore the confidence of a patient. Decisively going beyond products, we bring ongoing training and education that not only helps professionals stay relevant, it puts them ahead of the game.

We actively advance and connect dentists and lab technicians – the people who are the backbone of the industry.

Advance. Connect. Together.

www.3shape.com



Contact your reseller regarding availability of 3Shape products in your region.

*References/data available upon request from 3Shape

86660064 TRIOS Brochure 2024 EN Global

© 3Shape A/S, 2024. All rights reserved.
The 3Shape name and logo and/or other trademarks mentioned herein are trademarks of 3Shape A/S, registered in EU and other countries.

3shape 