



because clean matters.



## SweepZone® Ultrasonic Cleaning Systems



SweepZone 200 shown w/accessories

All SweepZone Ultrasonic Cleaning Systems operate with a special cleaning wave of +/- 2KHz energy that literally sweeps back and forth throughout the ultrasonic tank. The innovative L&R SweepZone technological system incorporates Pulse Width Modulation (PWM), which enables the machines to sense cleaning loads, solution levels and voltage fluctuations.

SweepZone Ultrasonic Cleaning Systems automatically adjust to different conditions in order to create more uniform cleaning power by eliminating hot spots and weak areas in the tank. This results in unsurpassed cleaning performance with every application.

User friendly digital touch pad controls make it easy and convenient to set any of the five preset timer settings from a range of six to sixty minutes. Durable stainless steel housing and tank are corrosion resistant and long lasting.

Product Code	Description	Input Power Average watts 117 VAC 60 Hz	Amps	Overall Dimensions* Inches Centimeters (L x W x H)	Shipping Weight Pounds Kilograms	Tank Capacity Gallons & Liters	Tank Internal Dimensions Inches Centimeters (L x W x D)
AG902	<b>SweepZone 200</b> w/Timer & Drain, Stainless Steel Cover	100	0.812	12.68 x 6.88 x 8.75 32.2 x 1.75 x 22.2	12 5.4	1.06 4.0	11.75 x 6 x 4 29.8 x 15.2 x 10.2
AG910	<b>SweepZone 310</b> w/Timer & Drain, Stainless Steel Cover	400	3.419	16.5 x 10 x 12.13 41.9 x 25.4 x 30.5	26 11.8	3.25 12.3	15.5 x 9 x 6 39.4 x 22.9 x 15.2
AG918	<b>SweepZone 1000</b> w/Timer & Drain, Stainless Steel Cover	500	4.274	17 x 19 x 20 43.2 x 48.3 x 50.8	101 45.8	11.62 44.0	14 x 16 x 10 35.6 x 40.6 x 25.4
AG919	<b>SweepZone 1000</b> w/Timer, Drain, Heat, Stainless Steel Cover	500	4.274	17 x 19 x 20 43.2 x 48.3 x 50.8	101 45.8	11.62 44.0	14 x 16 x 10 35.6 x 40.6 x 25.4

All systems operate at 43kHz – 45kHz sweep range with a sweep rate of either 3 or 83 Hz. All systems are available with international voltages

**All L&R SweepZone Ultrasonic Cleaning Systems are backed by an unprecedented 2 ½ year warranty on parts & labor.**



SweepZone 310



SweepZone 1000

Visit our website to see our full line of quality solutions and accessories

# SweepZone® Recessed Ultrasonic Cleaning Systems



L&R's Recessed Systems with digital touch pad controls easily install into existing countertops and work areas to save valuable workspace.

These models provide an efficient and convenient method for instrument cleaning & processing. SweepZone recessed systems utilize simple inline connections that facilitate installation with most cabinet and Steri-Center Work Stations.

(Cabinet Not Included)

Product Code	Description	Input Power	Output	Overall Dimensions*	Shipping Weight	Tank Capacity	Tank Internal Dimensions
		Average watts 117 VAC 60 Hz	Frequency	Inches Centimeters (L x W x H)	Pounds Kilograms	Quarts & Liters	Inches Centimeters (L x W x D)
AG 916	<b>SweepZone 310R</b> w/Remote Digital Timer, Drain, Mounting Rim	400	43 KHz +/-2KHz 83Hz	16.25 x 10.13 x 6.25 41.3 x 25.7 x 15.9	26 11.7	13 12.3	15.5 x 9 x 6 39.4 x 22.9 x 15.2
AG 915	<b>SweepZone 360R</b> w/Remote Digital Timer, Drain, Mounting Rim	300	43 KHz +/-2KHz 83Hz	16.25 x 10.13 x 6.25 41.3 x 25.7 x 15.9	22 9.9	14.5 13.6	11.5 x 9.25 x 8 29.2 x 23.5 x 20.3
AG 917	<b>SweepZone 650R</b> w/Remote Digital Timer, Drain, Mounting Rim	500	43 KHz +/-2KHz 83Hz	16.25 x 10.13 x 6.25 41.3 x 25.7 x 15.9	38 17.2	26 24.6	19.75 x 11.75 x 6.5 50.2 x 29.8 x 16.5



Product Code	Stainless Steel Solid Side Draining Baskets
18666	200 (Full Size)
13518	310, 310R (Full Size)
18202	310, 310R (1/2 Size)
18676	360, 360R (Full Size)
17954	650, 650R (Full Size)



S200 Accessory Kit (PC 659)



310 Cassette Rack (PC 26665)



Solid Side Basket (PC 13518)

Product Code	Instrument Cassettes & Racks
17752	Large Cassette 310, 310R, 650, 650R
19368	Cassette Rack for 360R, Holds 5 Small Cassettes
17583	Accessory Box
17582	Bur Box (Holds 32 Burs)
26665	310, 310R Cassette Rack (Holds 3 Large Cassettes or 6 Small Cassettes)
15185	650, 650R Cassette Rack (Holds 3 Large Cassettes or 6 Small Cassettes) 2 Racks Fit 1 Unit

Product Code	Accessory Kits
659	200 - 3/4 Size Solid Side Basket, Positioning Cover, 2-600 ML Beakers Complete, Bur Box, Bur Tray
535	310, 310R - 1/2 Size Solid Side Basket, 1/2 Size Positioning Cover, 2-600 ML Beakers Complete, Bur Box, Bur Tray
1214	360, 360R - Solid Side Basket, Positioning Cover, 4-600 ML Beakers Complete, Bur Box, Bur Tray
1132	650, 650R - 1/2 Size Mesh Basket, 1/2 Size Positioning Cover, 4-600 ML Beakers Complete, Bur Box, Bur Tray

Note: Recessed housings are vinyl coated steel.



L&R Manufacturing Company  
577 Elm Street, P.O. Box 607  
Kearny, NJ 07032-0607 USA  
(P) 866.567.8811 / 201.991.5330 (F) 201.991.5870  
www.LRultrasonics.com info@LRultrasonics.com