



500 XLS

SURGICAL CHAIRS & TABLES FOR EYE- & MICROSURGERY



E

ENJOY QUALITY AND DESIGN



500 XLS



For IVOM and ENT treatments

- Versatile design to support a wide range of application: OMF, dental, ophtalmmic surgery
- Outstanding price-performance ratio
- 4 powerful drives for adjusting the back, foot rest and seat elements as well as the overall height
- With equipotential bonding connection
- Mono lifting column designed for a maximum patient wieght of up to 200 kg
- Convenient push handles and smooth-moving castors
- Comfortable upholstery in anthracite standard shade
- Easy entry and exit
- Comfortable armrests, removable



Item No. : 1500 0005

OUR SPECIALIST



500 XLS



Entering / exiting

The chair's low entry height allows the patient to get easily seated in a comfortable position. After completion of the treatment, the treatment chair can be easily returned to its reset position by stepping the foot keypad. The patient is able to leave the chair safely and comfortably.



3 MIN



Surgery position

The patient is moved into the perfect position for any type of treatment.

Patient positioning

The memory settings allow rapid and safe patient positioning at the push of a button via hand or foot control.

FUNCTIONALITY IN A NUTSHELL



Precise positioning of the head part

Headrest

Multi-functional headrest element featuring two-dimensional and longitudinal adjustment. Horseshoe shape also allows treatments in prone position.



Braking system

Mechanically. Easy to handle. Rugged moving gear featuring 2 central braking steering castors and 2 rigid castors for high stability combined with easy maneuverability.



Mains-free operation

Power supply-independent battery box with integrated emergency shut-down switch. Operating time per battery charge approximately 100 OPs.



Easy handling




Positioning by hand control

Multifunctional hand control with 8 memory positions, auto-run function und reset key.

Positioning by foot control

Removable foot keypad for back, seat and height adjustment, including on-off button, 3 memory positions, and a reset key.



Transport



Anesthesia



Treatment / surgery 1



Flat position

DEVELOPED FOR **MAXIMUM STABILITY**

Also think about the little things - Upgrade your OR to perfection.



assistTrend



assistPro H



easyPro F



supportLine



UFSK-INTERNATIONAL OSYS GmbH
Kirchhoffstraße 1
93055 Regensburg
Germany

MADE IN GERMANY

Design | Construction | Manufacturing

Phone
Fax
Mail
Web

+49(0)941 78862 15
+49(0)941 78862 35
info@ufsk-osys.com
ufsk-osys.com