

Evanesce™ Flow

Flowable Restorative



Today's Minimally Invasive Dentistry Demands More from Flowable Composites.

In addition to its traditional use as a base liner under posterior composite restorations, the ability to provide a strong, durable, and esthetic conservative restoration using a flowable composite is an invaluable asset to any dental practice.

Evanesce™ Flow has been engineered to complement the **Evanesce** Universal Restorative composite line and features the very same chemistry; resulting in amazing esthetics, high strength and preferred handling. The shades and opacities available in **Evanesce Flow** mirror the applications that **Evanesce** is designed for, and provide an excellent shade and opacity match to **Evanesce** and the VITA Classic shade guide, enabling it to be used on its own or in combination with **Evanesce**. **Evanesce Flow** moves easily when manipulated without sticking to the instrument causing pull-back. **Evanesce Flow** is thixotropic for improved control, is easily dispensed from the syringe or a single dose, and will not slump on a vertical surface – making it ideal for use in Class V's and small Class II's. Regardless of its application, its relatively low rate of shrinkage (3.6% compared with up to 5%), and high filler content (62% by weight and 48% by volume), make it truly universal for use anywhere in the mouth.

Whether You're Choosing for Shade, Location or Opacity — There's an Evanesce Flow For That!

Evanesce Flow is ideal for superficial Class I and II's due to its esthetics, ability to be adapted to the tooth and high strength. Class V's benefit from its no-slump handling and range of opacities and shades. **Evanesce Flow** is invaluable as a liner base in posterior restorations to support the dentinal/adhesive seal, preventing post-op sensitivity due to its low contraction shrinkage. Radiographically, **Evanesce Flow** is easily distinguishable from tooth structure and neighboring composite with a radiopacity of 200% Al.

Available in the most popular shades in Enamel, Universal and Dentin opacities, **Evanesce Flow** will easily satisfy the majority of your flowable composite applications. Versatile **Evanesce Flow** is ideal for small esthetic restorations in both the young as well as the aging populations; and for blocking out undesirable tooth structure or layering restorations where there are varying degrees of cavity depth.

Esthetic and strong, easy to express and manipulate, available in the most popular shades and opacities, **Evanesce Flow** is the ideal solution to minimally-invasive restorations and a vital component to your large composite cases.



Evanesce Flow is clearly distinguishable from composite on a radiograph for ease of diagnosis.

EVANESCE FLOW Shades and Opacities

ENAMEL SHADES (80% Opacity)					
A1	A2	A3	B1		
UNIVERSAL SHADES (85% Opacity)					
A1	A2	A3	A4	B1	BL1
DENTIN SHADES (90% Opacity)					
		A2			

Available in 1.7g Syringes or 0.25g Single Dose

EVANESCE FLOW 1.7g SYRINGE

Contains: 1 x 1.7g syringe, mixing tips, instructions/SDS

A1 Universal (323517)
 A2 Universal (323518)
 A3 Universal (323519)
 A4 Universal (323520)
 B1 Universal (323521)
 BL1 Universal (323525)
 A1 Enamel (323501)
 A2 Enamel (323502)
 A3 Enamel (323503)
 B1 Enamel (323505)
 A2 Dentin (323511)

Unique unit-dose compule with 40° pivoting tip



EVANESCE FLOW 0.25g SINGLE DOSE

Contains: 20 x 0.25g single dose, instructions/SDS

A1 Universal (323556)
 A2 Universal (323557)
 A3 Universal (323558)
 A4 Universal (323550)
 B1 Universal (323559)
 BL1 Universal (323555)
 A1 Enamel (323531)
 A2 Enamel (323532)
 A3 Enamel (323533)
 B1 Enamel (323535)
 A2 Dentin (323546)