



NEOBURR®

The greatest invention since sliced bread

More blade contact with the tooth means fantastic cutting ability without clogging up the blades.

We don't have a breaking point

Blended Neck® technology on most NEOBURR carbides provide superior strength and durability. Being constructed with a single piece of material eliminates the weld that is often the point which other burs separate. No breaking point here.



Testing, testing 1-2-3

We tested it but others did too. In 3rd party testing our 557² carbide was shown to be up to 70% stronger than the competition.

Nobody likes a chatter box.

With its unique construction NEOBURR provides a fast, smooth, chatter-free experience. This is good for you and your hand-piece.

"Love your product. Price and packaging, quality of cut, all exceeds my expectations."

E. Rothfield DDS,
Hollywood, FL