

# A.S.A.P.<sup>®</sup>

## All Surface Access Polishers



### More Shine. Less Time.

Diamond-impregnated **A.S.A.P.<sup>®</sup> Polishers** provide all surface access polishing and are ideal for use on all composite surfaces, lithium disilicate, bisacryls, most ceramics, and zirconia.


- The A.S.A.P. Pre-polisher (purple) is embedded with diamond particles (44 microns) and gently diminishes small surface defects without affecting anatomy, and prepares the surface for a final high gloss shine.
- The A.S.A.P. Final High Shine Polisher (peach) is embedded with small, 3-6 micron sized diamond particles to provide a life-like polish in as little as 20 seconds.
- Autoclavable up to 30 times, A.S.A.P. Polishers can withstand wear while providing an efficient and consistent polish use after use.

### All Surface Access Polishing



### Proven Performance

A recent university-based study,\*\* comparing the gloss produced by A.S.A.P. Polishers (Clinician's Choice<sup>®</sup>) and Enhance<sup>®</sup>/PoGo<sup>®</sup> (Dentsply Sirona) found that A.S.A.P. Polishers produced a higher gloss results, after 30 seconds, than Enhance and Pogo did after 90 seconds.



\*\*Data on file

“A.S.A.P. Polishers save chair time by quickly polishing porcelain and composite restorations in **less than a minute.**”

- Danièle Larose, DMD, AAACD



Dentistry and photography courtesy of Danièle Larose, DMD, AAACD.



Pre-op



A.S.A.P. Pre-polisher



A.S.A.P. Final High Shine Polisher



After

**CLINICAL Instant Replay** | In a rush? Watch this quick technique tip! | POWERED BY Hands-OnLine LIVE

Scan to watch

**Achieving a Superior Polish with A.S.A.P. All Surface Access Polishers**

Bob Margeas, DDS

**CLINICAL Instant Replay** | In a rush? Watch this quick technique tip! | POWERED BY Hands-OnLine LIVE

Scan to watch

**Composite Finishing Technique for an Esthetic Result**

Bob Margeas, DDS



Request a complimentary sample:  
[sample.clinicianchoice.com/asap-polisher-family](http://sample.clinicianchoice.com/asap-polisher-family)



<p><b>A.S.A.P. STARTER KIT (130917)</b>                  Contains: 3 x Pre-polishers Small, 3 x Pre-polishers Large, 3 x Final High Shine Polishers Small, 3 x Final High Shine Polishers Large</p>	<p><b>A.S.A.P. TRIAL 2-PACK</b>                  Contains: 1 x Pre-polisher, 1 x Final High Shine Polisher                  Small (130910)                  Large (130911)</p>
<p><b>A.S.A.P. REFILL 6-PACK</b>                  Contains: 3 x Pre-polishers, 3 x Final High Shine Polishers                  Small (130913)                  Large (130912)</p>	<p><b>A.S.A.P. BULK PACK SMALL PRE-POLISHER (130927)</b>                  Contains: 12 x Pre-polishers Small</p>
<p><b>A.S.A.P. PRE-POLISHER 3-PACK</b>                  Small (131011)                  Large (131012)</p>	<p><b>A.S.A.P. BULK PACK LARGE PRE-POLISHER (130928)</b>                  Contains: 12 x Pre-polishers Large</p>
<p><b>A.S.A.P. FINAL HIGH SHINE POLISHER 3-PACK</b>                  Small (132013)                  Large (132014)</p>	<p><b>A.S.A.P. BULK PACK SMALL FINAL HIGH SHINE POLISHER (130925)</b>                  Contains: 12 x Final High Shine Polishers Small</p>
	<p><b>A.S.A.P. BULK PACK LARGE FINAL HIGH SHINE POLISHER (130926)</b>                  Contains: 12 x Final High Shine Polishers Large</p>

To view the latest news, education and more, visit:  
[www.clinicianschoice.com](http://www.clinicianschoice.com)

