

DIGITEST³

Digital Pulp Vitality Tester

Product reviewed by



Read Style Italiano's review here

Identify the Source of Your Patient's Pain

- ▶ Single tooth testing.
- ▶ Intuitive single button operation.
- ▶ Multiple test settings to accommodate your personal preference.
- ▶ Four autoclavable probes (included) promote easy access to all tooth surfaces.
- ▶ Standard 9V battery enables simple in-office replacement.



Dr. Giuseppe Chiodera, Style Italiano


Elevate Your Diagnostic Skills to the next Level

When percussion, ice, heat and radiographs are not enough, electric stimulation may be the diagnostic tool you need for determining the cause of the patient's discomfort. That's where the Digitest³ comes in.

How it works

The Digitest³ electrically stimulates the tooth to test for vitality, resulting in a value ranging from 0 to 64. Comparing the response to the stimulus of the suspected tooth with a known healthy tooth's response to an identical stimulus helps you:

- ✓ assess tooth vitality
- ✓ identify the source of your patient's pain
- ✓ test and follow up traumatized teeth.



SKU #	PRODUCT DESCRIPTION
D655	Digitest³ Pulp Vitality Tester (4) Autoclavable Probes; 1-standard long, 1-standard short, 1-labial, 1-lingual, (1) Ground Clip, (1) Lead Wire, (1) 9V Battery
D636	Precision Probe Set (2) Probes; one labial, one lingual
D625	Standard Metal Probe Set (2) Probes; one short & one long
D631GC	Pulp Tester Ground Clip & Lead Wire (2) Clips, (1) Wire



For information and orders contact infoeurope@parkell.com | +46 708 593481

